

Batch Certificate For Research Use Only

PRODUCT INFORMATION AND QUALITY CONTROL

NAME OF PRODUCT	SID-Plasma 2x (human) DNAfree		
DESCRIPTION	Human serum proteins in common plasma concentrations, electrolytes, EDTA		
CATALOG NUMBER	SID-0001		
BATCH NUMBER	SID-XXXX		
MANUFACTURING CONDITIONS	<ul style="list-style-type: none"> • Manufactured und sealed in class 2 safety cabinet • Bottled with qualified OPUS® Dispenser • At room temperature 		
PACKAGE SIZE AND TYPE	<ul style="list-style-type: none"> • 2D barcoded tube with screw cap • Material: Polypropylen (PP) 		
DATE OF MANUFACTURE	01 Oct 2018		
EXPIRY DATE	01 Oct 2019		
CONCENTRATION	2x		
NOMINAL VOLUME	2,5 ml		
QUALITY	Sterile, DNAfree		
STORAGE CONDITIONS	+ 2-8 °C		
MANUFACTURING AND QUALITY CONTROL SITES	SensID GmbH Schillingallee 68, 18057 Rostock, Germany		
TEST METHOD AND ACCEPTANCE CRITERIA	Quality Control	Test Method	Acceptance criteria
	Quality (DNAfree)	Real-Time PCR Investigator Quantiplex Pro Kit (Qiagen)	< 0,5 pg/μl ¹
	¹ LOD		

RESULTS OF ANALYSIS	Result	PASS/FAIL
	Quality	< 0,5 pg/ μ l
COMMENTS/REMARKS	Additional information that may be relevant to the importer and / or inspector. Compliance of the batch certificate (e.g. specific storage or transport conditions).	

Name and position/title of Person authorising the batch release:

Mr. Björn Nowack

Date of batch release: 02 Oct. 2018

Signature batch release: Björn Nowack